

LAPAROGUARD[®]

STERILE DISPOSABLE LAPAROSCOPIC INSUFFLATION FILTER



- Disposable Bacterial/Viral filter
- 0.027 micron filter
- Ultraclean made from styrolux
- 99.999% Filtration efficiency
- Prevents the risk of cross contamination

ORDERING INFORMATION

Ref No	Qty / Minipack
6010148	25



 **MORRISONS**
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com