



Lifecare[®]

SHIPPER :

**200 / Box
800 / Carton**

VISOR MASK

- 4-Ply Tie-On Mask with Anti Fog Optical Visor
- 175 Micron PET Optical Visor
- Black Fabric Strip over Noseclip deflect light glare
- Long and soft PU Foam Strip under nose clip for additional user comfort and reduce fogging. Extra long noseclip.
- Strong Tie Bands, 40 Gsm, prevents tie snapping
- 90 cm Horizontal Tie Bands accommodates all head sizes and hair lengths
- Malleable PVC coated aluminium nose-clip prevents air leakage
- High Bacteria Filtration* (4-Ply: BFE>99%)
- Low Breathing Resistance* (4-Ply: Delta P < 49,0 Pa)



MORRISONS
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com