

ENTERAL FEEDING BAG

WITH TRANSFER SET



DESCRIPTION

- Specially designed for Gastro-Intestinal feeding.
- Provided with transparent and flexible drip chamber.
- Rigid neck with closure cap & hanger for easy filling and hanging.
- An Ice pouch is provided on the rear side of the bag for preserving the food longer time.
- Smooth roller clamp for regulating the flow rate.
- Graduated bag with markings.
- Available in 500ml & 1000ml.
- Sterile, individually packed.

ORDERING INFORMATION

Ref No.	Capacity	Qty / Mini Box	Qty / Master Carton
7010120	500ml	15	60
7010121	1000ml	15	60



MORRISONS
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com