

MOLIDRAIN[®] - IR

Interventional radiology procedure drainage bag



DESCRIPTION

- Compactable with all standards of Luer locking drainage catheters and can remove fluids using gravity.
- Secure the connection to the catheter, maintaining a closed system to reduce the chance of leakage or infection.
- Super smooth, kink-resistant tubing ensures unrestricted flow.
- Easy visualisation for when drainage bags require emptying
- The T' type closure provides quick and easy emptying of the drainage bag.
- The anti-reflux valve prevents retrograde contamination.
- A specially designed moulded handle facilitates easy carrying and hanging of the bag.
- Sterile, individually packed.

ORDERING INFORMATION

Ref No.	Qty / Mini Box	Qty / Master Carton
6010155	40	160



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com