

GUEDAL AIRWAYS



Product Features :

- Design of Hard bite block facilitates passage of a suction catheter.
- Rigid bite block to ensure airway integrity and avoids biting of tongue.
- Smoothly rounded edges for improved patient comfort.
- Colour coded for instant identification of sizes.

ORDERING INFORMATION

Ref No	Size	Qty/Box	Qty/Case
3010143	00	300	1200
3010144	0	300	1200
3010145	1	250	1000
3010146	2	150	600
3010147	3	150	600
3010148	4	150	600

