


SPIROMETER




Product Features :


Compact, high-quality, tear-resistant plastic, 3-ball respiration measurement system
Solid and non-shaky construction; removable bottom for deep cleaning and hygiene maintenance.
Deep breathing exerciser for complete breathing training - 3 chambers for different inhalation rates from 600 cc / s to 1200 cc / s.
Helps achieve optimal lung capacity and restore interrupted breathing patterns.
Complete instructions supplied; view video in images next to the ball; colors may vary




01 Place your mouth around the mouthpiece and inhale



02 Inhale well to lift the balls to the top



03 Inhale for as long as you can comfortably



04 Remove the mouthpiece and breathe out properly. Repeat the process

