


SUCTION CATHETER



Product Features :

- Non-toxic, non-irritant medical grade PVC.
- Distal end is open with one lateral eye.
- Frosted surface tubing for smooth intubation.
- Colour coded for instant size identification.

ORDERING INFORMATION

Ref No	Size	Qty / Box	Qty / Case
3010101	5fg	100	1300
3010102	6fg	100	1300
3010103	8fg	100	1300
3010104	10fg	100	1300
3010105	12fg	100	1300
3010106	14fg	100	1300
3010107	16fg	100	1300
3010108	18fg	100	1300

