


QUEENCATH[®]

URETHRAL CATHETER



Product Features :

- Especially designed for short term bladder catheterization of females.
- Manufactured from non-toxic, non-irritant, soft medical grade PVC.
- Atraumatic, soft, rounded, closed tip with two lateral eyes for efficient drainage.
- Universal funnel shaped connector for safe connection to urine bags.
- Radio opaque line throughout the length for x-ray visualization.
- Colour coded plain connector for easy identification of size.
- Sterile individually packed in peelable pack.

ORDERING INFORMATION

Ref No	Size	Qty/Box	Qty/Case
5010110	8fr	500	2000
5010112	10fr	500	2000
5010113	12fr	500	2000
5010114	14fr	500	2000
5010115	16fr	500	2000
5010116	18fr	500	2000
5010117	20fr	500	2000

