

# Queencath CIC Kit

**50 Female Urinary Catheters with a CIC cloth Pouch**



- Used for short-term bladder catheterization.
- Radio opaque line provided throughout the catheter for X-ray visualization.
- Frozen surface tubing for super smooth intubation.
- Atraumatic, soft rounded, closed tip with two lateral eyes for efficient drainage.
- The universal connector at the proximal end for easy connection to the urine bags.
- Available Sizes (FR) : 8,10,12,14,16,18,20.
- Kit contains 50 catheters & a CIC cloth pouch.
- Sterile, individually packed.

